To the Heads of
the Regional Units / Local Companies

Name	Werner Schneider
Dept.	SP WH
Telephone	09131/84-4394
Fax	09131/84-8456
Mail	schneider.werner@siemens.com
our Ref.	
Date	2003-10-24
Complaint- number:	SP-C-03.0285

Safety note

Dear colleagues,

we herewith ask you to pass the enclosed letter to all customers concerned in your sales region.
This letter should be filed in the register "safety" of the respective operating manual.

ATTENTION:

Please distribute immediately, by registered letter with reply coupon, the letter "Important Safety Information" (see enclosure) to the customers concerned .

If necessary, please have the enclosed customer letter translated technically correct into the local language and copy it yourselves. It is not permitted to change the content of the letter.

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Potential for injury due to Compression Plate Damage

Dear Customer,

The safety and comfort of your patients are extremely important to us. As a precaution, we would like to advise you of the following:

The normal aging process of the plastic in the compression plates (cradle) can be accelerated by high mechanical stress during compression, as well as the use of unsuitable cleaning agents.

As explained in the operating instructions, compression plates are considered routine consumables subject to wear and tear.

A letter from a customer has shown that extreme damage in the cradle can result in patient injuries.

Our recommendations for safe handling of compression plates are as follows:

- Use the recommended cleaning agents only.
- Prior to use, check the cradle for any signs of damage in the area near the chest wall, for example, cracks.
The compression plate may not be used if damaged.
- Replace the compression plates if there are any signs of damage.

Please keep this safety notice with your operating instructions.

Sincerely,

Siemens Aktiengesellschaft

Gerhard Köbrunner

Dr. Klaus Wecker

Division
Special Systems

Head:
Holger Schmidt

Medical Solutions
Group Executive Management:
Prof. Dr. E. R. Reinhardt, Group
President
Dr. Hermann Requardt
Götz Steinhardt

Postal address:
Siemens AG
SP QP
P.O. box 3260
D-91050 Erlangen

Office address:
Allee am Röthelheimpark 2
D-91052 Erlangen
Telephone (09131) 84-0

Siemens Aktiengesellschaft · Chairman of the Supervisory Board: Karl-Hermann Baumann · Managing Board: Heinrich v. Pierer, Chairman, President and Chief Executive Officer · Members: Johannes Feldmayer, Thomas Ganswindt, Klaus Kleinfeld, Edward G. Krubasik, Rudi Lamprecht, Heinz-Joachim Neubürger, Jürgen Radomski, Erich R. Reinhardt, Uriel J. Sharef, Claus Weyrich, Klaus Wucherer
Registered Offices: Berlin and München · Commercial Registries: Berlin-Charlottenburg, HRB 12300; München, HRB 6684

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